



Request on Mass Spectra (Thermo Exactive Plus)

Name:.....Phone:.....Date:.....

Working group:.....

Sample code number:

			-						-		
--	--	--	---	--	--	--	--	--	---	--	--

Structure of expected product and reaction equation with all reactants and solvents:

Molecular weight: **Da** **Sum formula (Monomer):**

Mass range of interest: from to Da

weighed portion [mg] 0.5 mg 1.0 mg 2.0 mg mg (please mark)

Suitable Solvent(s): (solvents such as DMSO, DMF, pure water, and toluene are generally not suitable for MS measurements)

The sample is: air-sensitive moisture-sensitive light-sensitive

toxic harmless

Method of measurement:

ESI APCI ASAP LIFDI